

# Advantech AE Technical Share Document

<b>Date</b>	2021/7/30	<b>Release Note</b>	<input type="checkbox"/> Internal <input checked="" type="checkbox"/> External
<b>Category</b>	<input checked="" type="checkbox"/> FAQ <input type="checkbox"/> SOP	<b>Related OS</b>	
<b>Abstract</b>	What Camera Utility Should I Use for Quartz Camera(QCAM)?		
<b>Keyword</b>	Quartz Camera, QCAM, Camera Utility		
<b>Related Product</b>	QCAM, PCIE-1672, PCIE-1674, PCIE-1181, PCIE-1182		

## ■ **Problem Description:**


Users may need to use camera utility to adjust certain camera settings such as AOI, exposure time, gain, and so on from time to time to accommodate their vision applications. It is sometimes confusing to determine which camera utility to use as it will differ in accordance with the camera and the frame grabber used.

In this documentation, you will find information regarding which camera utility to use and where to get it step by step.

## ■ **Brief Solution - Step by Step:**

If you connect Quartz camera to PCIE-1672 or PCIE-1674:


1. Search “Pylon Viewer” on browser(or click on [here](#)) to download the latest camera software (e.g. pylon 6.2.0 Camera Software Suite Windows).
2. Run the installer and select a suitable user profile & camera interface to install.
3. Use pylon IP Configurator to set IP address and subnet mask for NIC and camera.
4. Use pylon Viewer to connect to camera, adjust camera setting, and capture live images.



5. (Optional, apply only when necessary)

*[Note]: If you find the image capture or parameter setting is somehow irresponsive or corrupted, continue to apply the following step 5 and install a special pylon license.*

● **Issue screenshot:**



Download the license pack from Advantech Quartz Camera Download (QCAM Optional License). Follow the instructions in the readme document to finish installing the license files.

Follow the step

1. Unzip the "License.zip" file.
2. For Windows copy the "\*.lic" and the "pylon.xml" files into the pylon Runtime installation folder where the "PylonBase\_XXX.dll" resides, e.g.:
  - for 64bit:
   
C:\Program Files\Basler\pylon X\Runtime\x64
  - for 32bit:
   
C:\Program Files\Basler\pylon X\Runtime\Win32
3. For Linux copy the "\*.lic" and the "pylon.xml" files into the pylon Lib installation folder where the "libpylonbase-XXX.so" resides, e.g.:
  - for 64bit:
   
opt\pylonX\lib64
  - for 32bit:
   
opt\pylonX\lib

Make sure it can capture live images.

